# Evolutis

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**JUST  
UNIC®**

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**Jack-Up** Shoulder implant  
for **Trauma**

# Jack-Up Shoulder implant for Trauma

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The JUST UNIC humeral implant is indicated for the repair of 3 and 4 fragments fractures of the proximal humerus independantly from the quality of the bone ((1) Report on 26 cases of Bilboquet device at mean FU of 25 months: "the peripheral support on the humeral head optimally counteracts any tendency toward displacement into varus".)

This device allows to realize in a fast and guided way, an anatomical and stable osteosynthesis of the complex fractures of the superior extremity of the humerus ((2) Prospective study of 22 Bilboquet cases at mean FU of 34 months: "Union was obtained in all fractures. There was no secondary tilting of the head, and no migration or pseudarthrosis of the tuberosities").

The implementation of the staple in the cephalic spongious bone provides a solid support onto which an effective distraction can be made. The sleeve and stem set, implanted in the humeral diaphysis allows to realize an adjustable distraction of the fracture site and to realize an anatomical restoration of the humeral metaphyso-cephalic arch. The reintegration of the tuberosities by lacing-up can then be made in a easy way.

This system, designed by the Professor L. Doursounian (St. Antoine teaching hospital, Paris, France) in the 1990s, allows to obtain a consolidation of the tuberosities and of the humeral head in a reproductive way. In cases of asymptomatic avascular humeral head necrosis, it is possible to convert the osteosynthesis into a prosthesis by changing the staple by a prosthetic head on the morse taper of the stem.

When the type of fracture justifies straightaway an indication of prosthesis (in particular in cases of fracture-dislocation) this device allows an easy adjustment of the height of the prosthetic head and insures an optimized repair of the tuberosities.

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(1) A new internal fixation technique for fractures of the proximal humerus - the Bilboquet device: a report on 26 cases.

L. Doursounian et Al. Dept of Orthopaedic Surgery, HEGP, Paris, France. J Shoulder Elbow Surg, July/August 2000: 279-288

(2) Complex proximal humeral fractures: A prospective study of 22 cases treated using the "Bilboquet" device, L. Doursounian, A.

Kilinc, B. Cherrier, G. Nourissat, Department of Orthopaedic Surgery and Traumatology, Saint-Antoine University Hospital, Paris, France, Orthopaedics & Traumatology: Surgery & Research (2011) 97, 58-66.

- + Solid osteosynthesis
- + Anatomic reduction of the fracture
- + Preservation of the humeral head
- + Consolidation of the tuberosities
- + Accurate positioning of the height of the humeral head
- + Efficient re-attachment of the tuberosities
- + Auto-graft in contact with the tuberosities (when used with the prosthetic head)

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#### Humeral head staple

Platform shaped to provide quality support to the femoral head

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#### Locking screw

Locks the stem to the sleeve at the proper height

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#### Humeral stem

Slides and docks into 1 of 4 positions to the humeral sleeve.

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#### Replacement humeral head

Option when the humeral head has lost its vascularisation. Hollow design for auto-grafting

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Jack-up forceps: adjust stem level to restore the humeral head height

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Designed and  
Manufactured in  
France

**CE** Evolutis  
Avenue de la Libération, 42720 Briennon, France  
0499 Tel : +33. (0)477.60.79.99 – Fax : +33. (0)477.60.79.90

[www.evolutisfrance.com](http://www.evolutisfrance.com)