Evolutis

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EBAUCHE

UNIC[®] STEMLESS

Surgical technique

Table of contents

| | |
|---|----|
| Disclaimer | 2 |
| Indications | 2 |
| Patient positioning | 3 |
| Surgical approaches | 3 |
| Operative steps | 4 |
| • Positioning of the humeral jig | 4 |
| • Introduction of the central pin for fixation of the resection guide | 5 |
| • Resection of the humeral cap | 6 |
| • Size selection of the humeral head | 7 |
| • Preparation of the humeral epiphysis | 8 |
| • Sizing of the STEMLESS component | 8 |
| • Preparation of the humeral anchorage | 9 |
| • Trail and reduction with the trial implants | 9 |
| • Implantation of the final STEMLESS implant | 10 |
| • Implantation of the final humeral head | 10 |
| Snapshot of the instrument set | 11 |
| References | 12 |

EBAUCHE

Disclaimer

This document is intended to be read only by experienced orthopaedic surgeons familiar with the surgical implantation of shoulder arthroplasty, and by individuals related to or acknowledged by the Evolutis company.

This publication is intended as the recommended procedure for using the Evolutis UNIC STEMLESS shoulder implants. It offers guidance only. Evolutis is the manufacturer of the device. As such and claiming no medical skill, Evolutis does not recommend a specific use of a product or a technique, therefore each surgeon should consider the particular needs of the patient and make appropriate adjustments where necessary. For any additional information related to the products, the indications and contra indications, the warnings and precautions of use, and the adverse effects, please refer to the INSTRUCTION FOR USE leaflet included in the packaging of implants. For further advice please contact your local representative.

Brand and products are protected for copyrights. It is strictly forbidden unless authorized by the Evolutis company, to duplicate or copy whole or part of this document.

Indications

The UNIC STEMLESS implantable device is designed to treat centered omarthrosis pathologies at the shoulder joint of physically young and active patients.

The UNIC STEMLESS implantable device is specially designed for fixation in the cancellous bone of the proximal humerus. It can be used as an hemi-arthroplasty of the shoulder joint, or associated with a glenoid anatomic resurfacing with the anatomic glenoid components of the conventional UNIC system.

The UNIC STEMLESS humeral component is designed for cancellous fixation without cement. The humeral component is made out of additive printing technology and results in nano-structured porous and interconnected trabeculae intended to favour in-depth osseointegration. The additional proximal HA coating enhances the secondary fixation process.

Should the UNIC STEMLESS be revised, the implant has been designed to facilitate the introduction of bone chisels along the 4 intra-osseous flanges, leaving the humeral bone near-genuine for the implantation of a conventional stemmed primary humeral implant.

Patient positioning

The patient should be positioned in a half upright position

The body of the patient on the table must allow for the operated arm to be free of the table edge and to be manipulated freely in extension and adduction without hindrance

Ideally the whole shoulder should be free and not hindered

A lateral support should be installed to support the arm alongside the body

EBAUCHE

Surgical approaches


The two most common approaches for prosthetic shoulder surgery are the Delto-pectoral and the deltoid split approaches. Both approaches present advantages, the main difference between the two approaches is the exposure of the Glenoid.

- deltoid split

- + simplicity
- + easier exposure to the glenoid
- + tuberosity repair is easier
- + sub-scapularis preserving approach
- cut through the deltoid muscle
- if the cuff is intact, the exposure is compromised

- delto-pectoral

- + familiarity, hence most common approach
- + respect of the deltoid muscle and of the supra-spinatus
- + exposure of the axillary nerve possible
- + humeral preparation
- subscapularis cut & repaired: increased risk of subscap rupture
- exposure of the glenoid may be compromised

a : delto-pectoral approach
b : deltoid split (McKenzie) approach

The preference between delto-pectoral and deltoid split approach is usually related to the necessity of properly exposing the glenoid bone.

In general the delto-pectoral approach will be preferred for a total anatomic shoulder prosthesis, while the deltoid split will be preferred for a total reverse shoulder prosthesis.

Positioning of the humeral jig

The instruments used in this chapter are

Humeral jig

E38 006

Check for positioning of both anterior and posterior teeth at the margins of the cervical neck (calotte céphalique). Lock the humeral jig in this position. Turn humeral jig perpendicular to the cervical neck line and read resection values on the jig. Check that resection height is consistent with the A/P value.

EBAUCHE

Introduction of central pin for fixation of the resection guide


The instruments used in this chapter are

AO quick connect for pin E38 037
Ø2.5 L.70mm pin E38 009

Place the AO quick connect for pin on the power tool.
Introduce one pin (Ø2.5 L.70mm) on the AO quick connect.
Introduce the pin into the humeral jig.
Press the lateral knob on the humeral jig, and pull the humeral jig out of the pin.

EBAUCHE

Resection of the humeral cap

The instruments used in this chapter are

| | |
|----------------------------|---------|
| Humeral resection guide | E38 008 |
| Resection level controller | E38 018 |
| AO quick connect for pin | E38 037 |
| Ø2.5 L.70mm pin | E38 009 |

Introduce the humeral resection guide on the center pin in the center hole marked "0".

Hold manually the humeral guide parallel to the humeral neck line.

Introduce the resection level controller into the slot of the resection guide and check for the medial and lateral levels of resection.

If required, modify the orientation of the resection guide.


Once the guide is in the correct "anatomic" position, introduce a second Ø2.5 L.70mm pin into the more medial "0" hole of the resection guide and drill into the humerus.

Lock the position of the resection guide with a third Ø2.5 L.70mm pin introduced in the more lateral and convergent hole of the resection guide.

Resect the humeral cap.

EBAUCHE

The instruments used in this chapter are
 Vernier caliper E38 036

Control the thickness of the resected humeral cap with the Vernier caliper. For the real value of the resection, use the sum of value read on the Vernier caliper plus the thickness of the saw blade (1.3mm). In illustration x the real resection value is: 13.3 (caliper) + 1.3 (saw blade) = 14.6mm

According to the measures read on the humeral jig (page y), the humeral head will be less than 46mm. The minimal corresponding thickness is therefore between 15mm (head size 44) and 17mm (head size 46).

In the case the head size 46 is confirmed at the following steps, a correction of the resection level may be necessary to approximate a thickness value of 17mm. To achieve an additional 2mm resection, remove the lateral convergent pin, re-position the resection guide on the 2 remaining pins and in the holes of the resection guide marked "2". Resect the excess of bone.

Size selection of the humeral head

EBAUCHE

The instruments used in this chapter are
 Humeral head trials E38 020 to E38 032

Select the humeral head trials corresponding to the initial measures (in our example 44 and 46), and position the humeral head trial on the resected humeral cap and/or on the resected humeral epiphysis. The humeral epiphysis will usually show a larger size than the humeral cap. Taking into consideration that the humeral cap may not be spherical, the surgeon needs to assess the best overall size adaptation, but a slightly smaller head is to be preferred to a slightly larger head.

In our example, the size 46.5 (blue trial) appears to have a better adaptation to both the cap and the epiphysis than the size 44 (pink trial). The larger size 49 (orange trial) is obviously too large and overhangs around the cap too much.

Preparation of the humeral epiphysis

EBAUCHE

The instruments used in this chapter are


| | |
|--------------------------|--------------------|
| Humeral head trials | E38 020 to E38 032 |
| Head centering sleeve | E38 019 |
| AO quick connect for pin | E38 037 |
| Ø3.5 L.70mm pin | E38 009 |

Introduce a head centering sleeve into the selected humeral head trial.
Place the humeral head trial at best on the humeral epiphysis.
Adjust the AO quick connect for pin on the power tool.
Introduce one pin (Ø3.5 L.70mm) on the AO quick connect.
Drill into the centering sleeve until contact with the lateral cortex of the humerus.
Remove the centering sleeve and the humeral head trial.

Sizing of the STEMLESS component

The instruments used in this chapter are

| | |
|------------------|--------------------|
| Centering plates | E38 011 to E38 015 |
|------------------|--------------------|

Select a centering plate of diameter 10 to 12mm less than the selected humeral head.
Introduce the centering plate on the Ø3.5 L.70mm pin.
The centering plate is crenelled to show 2 dimensions:
- the inner diameter (bottom of crenels) of the centering plate shows the overall diameter of the final implant. This dimension should remain at distance of the cortical bone.
- the outer diameter (top of crenels) of the centering plate provides an information regarding the necessary distance between the implant and the cortex : if one of the top sides of the crenels overhangs outside the humeral cortex, than the size is too large.

Preparation of the humeral anchorage

The instruments used in this chapter are

| | |
|-----------------------------|--------------------|
| Cannulated drill | E38 017 |
| Conformators Ø25 to Ø38 | E38 001 to E38 005 |
| Impaction handle with anvil | E38 033 |
| Modular trial taper | E38 007 |

Adapt the cannulated drill to the power tool (small AO connect.) and drill the center imprint up to the hilt of the reamer. Remove the drill.


Select the conformator of the sizes selected at the previous sizing of the component step.

Assemble the conformator with the impaction handle.

Introduce the conformator on the guiding pin and check for orientation of the flanges: orientate the flanges (shown with the "X" marks on the conformator) away from the most fragile regions of the cortex such as the quadricipital groove. In general, an orthogonal positioning of the flanges (in the A/P and M/L orientations) will be acceptable.

Impact firmly until the conformator comes in contact with the humeral cut.

EBAUCHE

Trial and reduction with trial implants


The instruments used in this chapter are

| | |
|---------------------|--------------------|
| Modular trial taper | E38 007 |
| Humeral head trials | E38 020 to E38 032 |

Remove the impaction handle and the centering pin.

Note: if total shoulder arthroplasty, leave the conformator on the humerus and proceed to the glenoid preparation.

Adapt the modular trial taper on the conformator. Position the trial humeral head on the taper. Reduce the shoulder joint and assess the mobility and stability. In case of joint laxity, from humeral head of size 44, 2 thicknesses are available.

Implantation of the final implant

The instruments used in this chapter are
 Impaction handle with anvil E38 033

Impaction handle with anvil E38 033
 Adapt the impaction handle on the trial taper and remove both the taper and the conformator.
 Ask for the final STEMLESS implant to be given and thoroughly check the size before opening.
 Open the sterile pack and leave the STEMLESS implant in the packaging foam.
 Screw the impaction handle directly on the STEMLESS implant.
 Position the implant on the humeral epiphysis and check for the proper orientation of the flanges.
 Impact the STEMLESS implant until the proximal ring of the implant comes in contact with the cut.
At this stage, do not impact the ring below the level of the cut.

EBAUCHE


Implantation of the final humeral head

The instruments used in this chapter are
 Impaction handle with anvil E38 033
 Head pusher tip E38 034

Ask for the final humeral head implant to be given and thoroughly check the size before opening.
 Open the sterile pack and seize the implant by hand.
 Position the implant directly on the taper of the STEMLESS implant.
 Assemble the head pusher tip on the impaction handle.

Impact the head on the taper.
 Impact firmly until the lower edge of the head comes in contact with the humeral edge.

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